


## DSP Crosslinker Protocol and Product Information Sheet

Product Category:	Homobifunctional Crosslinkers
Catalog Number(s):	<a href="#">c1106-100mg</a> , <a href="#">c1106-1gm</a> , c1106-custom
Product Name:	DSP Crosslinker
Alternative Name(s):	DSP; Lomant's Reagent; 3,3'-Dithiodipropionic acid di(N-hydroxysuccinimide ester); Di(N-succinimidyl) 3,3'-dithiodipropionate; Dithiobis-succinimidyl propionate
CAS Number:	57757-57-0
Chemical Formula:	C <sub>14</sub> H <sub>16</sub> N <sub>2</sub> O <sub>8</sub> S <sub>2</sub>
Molecular Weight:	404.42
Spacer Arm Length:	12.0 Å
Storage:	Upon receipt store at -20°C (shipped at ambient temperature). Protect from moisture (i.e. humidity); blanket under desiccated, inert gas.



### General DSP Crosslinking Protocol

1. Allow vial of DSP Crosslinker to fully equilibrate to ambient temperature before opening to prevent condensation inside the vial (DSP is moisture-sensitive and will hydrolyze).
2. Prepare a 50 mM solution of DSP crosslinker, by dissolving 10mg DSP in 495 µL of dry DMSO or dry DMF solvent. (Note: The NHS esters of DSP are susceptible to hydrolysis. If aliquoted into an aqueous solution and stored in it, then the DSP will degrade.)
3. Using a 20-fold molar excess approach (20:1 Crosslinker:Protein), add crosslinker solution to the protein sample in non-amine containing buffer (i.e. 25 mM Sodium Phosphate, pH 7.4), so that the final crosslinker concentration is between 0.5 to 5 mM. Optimal pH range is from 7 to 9.
4. Allow the sample to react at room temperature for 30-45 minutes. Allow slightly longer if sample must be kept on ice (recommended 2-3 hours). This reaction rate is not highly temperature sensitive.
5. Quench any unreacted DSP with 25 mM to 200 mM Tris, pH 7.4. Allow to react for 10-15 minutes at room temperature.
6. Desalt sample to remove unreacted DSP crosslinker (i.e. gel filtration, dialysis, etc.).

### Intracellular DSP Crosslinking Protocol

1. Remove media by washing cells twice with non-amine containing buffer (i.e. 25 mM Sodium Phosphate, pH 7.4).
2. The crosslinking solution as noted in steps 1 and 2 above, then add the crosslinker solution to the cells in a final concentration of ~2 mM.
3. Incubate the reaction mixture for 30-45 minutes at room temperature or for 2-3 hours if on ice.
4. Quench any unreacted DSP with 25 mM to 200 mM Tris, pH 7.4. Allow quenching reaction to proceed for 10-15 minutes at room temperature.

### Reference:

Wong, S.S. (1993) CRC Chemistry of Protein Conjugation and Crosslinking. CRC Press, Boca Raton, Florida.

ProteoChem \* 5570 W. 290 N. Suite 2 \* Hurricane, UT 84737  
Toll Free: 888.501.2436 \* Phone: +1.720.239.1630 \* Fax: +1.720.239.1631  
Technical Support E-mail: [ts@proteochem.com](mailto:ts@proteochem.com)  
[www.proteochem.com](http://www.proteochem.com)